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#### Cergy - France

+33 (0)1 34 20 78 00 tel.  
+33 (0)1 34 20 78 78 fax  
bioseprainfo@pall.com

Filtration. Separation. Solution.<sup>SM</sup>

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Rev. 02/2006

# Ultroser<sup>®</sup> HY and Ultroser SF

## Serum Substitutes for Cell Culture

### ■ Presentation

**Ultroser HY** and **SF** are serum substitutes successfully used to replace FBS (Fetal Bovine Serum) in cell culture media.

**Ultroser HY** serum substitute is used for growth of hybridoma and lymphoblastoid cells. It is available as a lyophilized powder corresponding to 20 mL of substitute, and in frozen form, in 400 mL vials.

**Ultroser SF** serum substitute is used in the *in vitro* culture of many anchorage-dependent cells which require a low steroid content medium. It is available in frozen form, in 10 mL vials.

Their biological activity is 5 times that of FBS.

### ■ Preparation

#### 1 - Lyophilized product (HY): Solubilization

The lyophilized powder is solubilized in 20 mL sterile distilled pyrogen-free water (water for injection) as follows:

- Open the vial by breaking the tip of the capsule
- Introduce the appropriate volume of water either using:
  - a sterile pyrogen-free pipette after first having removed the stopper,
  - or a sterile pyrogen-free syringe by injection through the stopper.

These two operations must be performed under sterile conditions.

Allow the lyophilisate to swell for 10 to 15 minutes. Solubilization is rapidly obtained by gently shaking the flask or by successive aspirations using a pipette or syringe.

The reconstituted solution should be clear and pink-coloured. For greater safety, the solution may be filtered through a sterilizing membrane of 0.22 µm pore diameter. This will not result in any loss of activity.

#### 2 - Frozen product (HY and SF): Thawing

Before use, leave the frozen vial at room temperature for a few hours. To speed up the thawing, cool water (30 °C) can be run over the vial continuously.

The reconstituted solution should be clear and pink-coloured. For greater safety, the solution may be filtered through a sterilizing membrane of 0.22 µm pore diameter. This will not result in any loss of activity. The product must be used up immediately after thawing.

### ■ Concentration

After reconstitution, **Ultroser G** and **HY** serum substitutes have a biological activity 5 times greater than that of FBS. They should be used at concentrations between 0.5 % and 4 % (2.5 % to 20 % of FBS). A concentration of 2 % (10 % of FBS) is generally used. To obtain a standard concentration of 2 %, mix 2 % **Ultroser** serum substitute with 98 % base culture medium. For each cell type, various concentrations from 0.5 % to 4 % will have to be evaluated to determine the optimal growth conditions. The effects on cell growth are not always proportional to the **Ultroser** concentration, nor in correlation with the results obtained with FBS in the same conditions of use.

### ■ Method of culture

#### 1 - Ultroser HY

##### Culture medium

Depending on cell requirements, different basal media such as Iscove (IMDM), RPMI, HAM's F12 or a HAM's F12 / RPMI mix may be used.

##### Cell transfer

Two techniques may be employed for sub-culture from a medium containing FBS:

- 1/ Direct transfer of the cells in the medium supplemented with **Ultroser HY** serum substitute,
- 2/ Or after an adaptation period by progressively reducing the quantity of serum, over a few passages:

- 1<sup>st</sup> passage: 50 % medium supplemented with serum (MSS) + 50 % medium supplemented with **Ultroser HY** serum substitute (MSU) in the optimum concentration.

- 2<sup>nd</sup> passage: 25 % MSS + 75 % MSU
- 3<sup>rd</sup> passage: 10 % MSS + 90 % MSU
- 4<sup>th</sup> passage: 100 % MSU

This second technique allows better adaptation of the cells to the new culture medium.

Before seeding the cells, check the pH of the final medium and, if necessary, adjust to  $7.2 \pm 0.1$ . Sodium bicarbonate or a synthetic buffer (HEPES) may be added to improve the buffering power of the medium.

## 2 - Ultrosor SF

### Culture medium

The culture medium may be the one generally used for culture with FBS (MEM, etc.). However, the use of an enriched medium such as HAM's F12, DMEM, IMDM, or RPMI is recommended.

### Culture support

The culture support is the same as that currently used with FBS. **Ultrosor SF** serum substitute contains adhesion factors. Therefore, it can be used for pre-incubation of the plastic dishes in the medium prior to adding the cells. However, for particular cell types which are difficult to grow, it may be necessary, to improve adhesion and differentiation, to coat the plastic dishes with one of the following products:

- poly-D-lysine,
- collagen,
- an extracellular matrix,
- adhesion factors such as fibronectin.

### Trypsinization

The trypsinization techniques for tissues and mono-layered cells are the same as those used with FBS. However, the trypsinization period should be as short as possible, to minimize cell damage. To remove any residual trypsin, it may be necessary to add a trypsin inhibitor such as soya bean extract, aprotinin..., at a

 **Ultrosor HY and SF serum substitutes contain a human protein that has been tested and found to be non-reactive to:**

- Hepatitis B surface antigens (HBsAg),
- HIV I, HIV II, anti-HCV antibodies,
- syphilis,

**in accordance with FDA-approved tests.**

**However, these products should be considered as potentially infectious agents during handling.**

concentration of 0.5 to 1 mg/mL, before dispersing the cells. If the cells still remain aggregated to each other, the cell suspension may be passed through a 40 µm cell sieve.

### Cell transfer

The technique for cell transfer is the same as that used for **Ultrosor HY** serum substitute (see Page 1).

## Storage

The lyophilized product (**Ultrosor HY**) is shipped at ambient temperature and must then be stored at 2 – 8 °C.

The product in frozen form (**Ultrosor HY** and **SF**) is shipped in dry ice and must be stored at -20 °C on arrival.

Once reconstituted, the solutions must be used up quickly. Otherwise they can be aliquoted and stored frozen at -20 °C for several weeks.

Frequent freezing and thawing must be avoided.

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© Pall Corporation 2006