INNOVATIVE TOOLS FOR

cytogenetics

*Synchro
set*

EK AMTS008

INSTRUCTION
FOR USE

CE

IVD

Euro  **clone**[®]

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SYNCHROSET: INSTRUCTION FOR USE

The "Synchroset" reagents perform very well with any lymphocyte culture protocol commonly used; anyway, much better results can be expected when they are used in combination with Euroclone's "Chromosome Kit" test tube.

These instructions are specifically designed for "Chromosome kit", but they will apply also to different procedures.

Lymphocytes culture

1. Aseptically add 0.5 ml of whole blood to the tube; tighten the cap and mix
2. Incubate at $37 \pm 2^{\circ} \text{C}$ (CO_2 not required)
3. After 48-56 hours, add 0,1 ml of Solution "A"; incubate overnight
4. The next morning, simply add 0,1 ml of Solution "B" and incubate for further five hours (it is not necessary to wash the cells)
5. Add 1-2 drops of Colchicin and incubate for 90 minutes*

* Should this incubation time result in a synchronization of cell cycle in a late phase of mitosis (late metaphase and/or anaphase), then we recommend to decrease the incubation time to 30, 45 or 60 minutes in order to determine the optimal period which will provide prometaphasic chromosomes.

PREPARATION OF SLIDES

Slides for chromosome observation can be now conveniently prepared following the Chromosome kit "fast protocol".

Recommendations

Unused aliquots of both Solutions "A" and "B" must be stored at $\leq -18^{\circ} \text{C}$, in the dark; repeated freezing and thawing will not alter the functionality of these reagents.

Chromosome which have undergone synchronization may need longer trypsinization times in order to obtain proper GTG or RHG banding.

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Euroclone S.p.A. - Life Sciences Division