

InVitroCare® SUPPLEMENTED IVF MEDIA

HTF with HSA Fertilization Medium

Formulated with EDTA, a stabilized form of glutamine, and 10% (v/v) HSA
180 day shelf life at 2-8°C
Cat # 2001-HSA - 20ml

IVC-TWO™ with HSA Cleavage Medium

Low Glucose and Phosphate free
Formulated with EDTA, Citrate, Glutathione, a stabilized form of glutamine, and 10% (v/v) HSA
180 day shelf life at 2-8°C
Cat # 2008-HSA - 20ml

IVC-THREE™ with HSA Blastocyst Medium

Phosphate free and elevated levels of Glucose
Supplemented with amino acids and vitamins and 10% (v/v) HSA
Formulated with Glutathione
180 day shelf life at 2-8°C
Cat # 2007-HSA - 20ml

InVitroCare Inc
8415 Progress Drive, Suite F, Frederick, MD 21701
800 786 1986 ~ 301 360 9660 ~ Fax 301 360 9661