

# InVitroCare® SPERM PREPARATION MEDIA

## SPERM WASH MEDIUM

*Contains 0.5% Human Serum Albumin*

HEPES-buffered formulation

18 month shelf life at 2-8°C

Cat # 2003 - 12 x 12ml or Cat # 2005 - 100ml or Cat # 2005-5 - 500ml

## SPERM CARE™

### **Ready-to-use gradient system for sperm preparation**

*Silane-coated colloidal silica particles suspended in HEPES-buffered HTF*

HEPES-buffered formulation

2 year shelf life at room temperature ~ After opening, store at 2-8°C

**UPPER LAYER (45%)**

Cat # 2221-50 - 50ml or Cat # 2221-100 - 100ml

**LOWER LAYER (90%)**

Cat # 2222-50 - 50ml or Cat # 2222-100 - 100ml

**100%**

Cat # 2220 - 100ml

**InVitroCare Inc**

8415 Progress Drive, Suite F, Frederick, MD 21701

800 786 1986 ~ 301 360 9660 ~ Fax 301 360 9661